

Pempidine

Catalog No: #S0296


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pempidine
Brief Description	Inhibitors
Purification	97.00%
Target Name	AChR antagonist; Others
Calculated MW	155.28
Formulation	C10H21N
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CN1C(C)(C)CCCC1(C)C
Pathways: Neuroscience; Others
Receptor: Nicotinic Receptor; ganglion
Boiling pt: 187-188°C
Melting pt: 187-188°C
Solubility: DMSO: 10 mM
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several months.

References

1. A. Spinks et al. Br. J. Pharmacol. Chem. Others. 1958, 13: 501-520.

Note: This product is for in vitro research use only