

## Lactitol monohydrate

Catalog No: #S0299


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Lactitol monohydrate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	362.33
Formulation	C <sub>12</sub> H <sub>24</sub> O <sub>11</sub> H <sub>2</sub> O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism

systemSMILES: C([C@@H]1[C@@H]([C@@H]([C@H]([C@@H](O1)O[C@H]([C@@H](CO)O[C@@H]([C@H](CO)O)O)O)O)O)OPathways: Other  
 sReceptor: OthersBoiling pt: 788.5CMelting pt: 95-98CSolubility: DMSO: 10 mMAppearance: White CrystalRemark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Kanters JA, et al. Acta Crystallogr C. 1990 Dec 15;46 ( Pt 12):2408-11.

Note: This product is for in vitro research use only