

Pargyline


Catalog No: #S0300

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pargyline
Brief Description	Inhibitors
Purification	97.00%
Target Name	MAO inhibitor
Calculated MW	159.23
Formulation	C11H13N
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CN(CC#C)CC1=CC=CC=C1Pathways:MetabolismReceptor:MAOBoling pt:228.4CMelting pt:Solubility:Water: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Patkar AA, et al. CNS Spectr. 2006, 11(5):363-75.

Note: This product is for in vitro research use only