

## Ethisterone

Catalog No: #S0301


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ethisterone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	312.45
Formulation	C <sub>21</sub> H <sub>28</sub> O <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine

systemSMILES:C[C@]12CCC(=O)C=C1CC[C@@H]3[C@@H]2CC[C@]4([C@H]3CC[C@]4(C#C)O)C

Pathways:Endocrinology/HormonesReceptor:ERBoiling pt:Melting pt:Solubility:DMSO: 0.05 mg/mL (0.16 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Herbai G, et al. Acta Pharmacol Toxicol (Copenh), 1971, 29(23), 177-193.

Note: This product is for in vitro research use only