

## Undecenoic acid


Catalog No: #S0306

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Undecenoic acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	184.28
Formulation	C11H20O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Metabolism systemSMILES:OC(=O)CCCCCCCC=CPathways:OthersReceptor:OthersBoiling pt:137°CMelting pt:23-25°CSolubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. EREAUX LP, CRAIG GE. Can Med Assoc J. 1949 Oct;61(4):361-4.

Note: This product is for in vitro research use only