


Description

Product Name	Actarit
Brief Description	Inhibitors
Purification	99.76%
Target Name	Others
Calculated MW	193.2
Formulation	C10H11NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:CC(=O)NC1=CC=C(CC(O)=O)C=C1Pathways:OthersReceptor:OthersBoiling pt:449.1?Melting pt:173-175CSolubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Canudas N, et al. Arzneimittelforschung. 2008;58(4):182-7.

Note: This product is for in vitro research use only