


Description

Product Name	Actarit
Brief Description	Inhibitors
Purification	99.76%
Target Name	Others
Calculated MW	193.2
Formulation	C10H11NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: CC(=O)NC1=CC=C(CC(O)=O)C=C1
Pathways: Others
Receptor: Others
Boiling pt: 449.1
Melting pt: 173-175°C
Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Canudas N, et al. *Arzneimittelforschung*. 2008;58(4):182-7.

Note: This product is for in vitro research use only