

Diathesin


Catalog No: #S0311

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Diathesin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	124.14
Formulation	C7H8O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OCC1=CC=CC=C1O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: 83-85°C
Solubility: Soluble in Water
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Michalski C, et al. J Biol Chem. 2008 Jul 11;283(28):19219-28.

Note: This product is for in vitro research use only