

Blutene chloride

Catalog No: #S0312


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Blutene chloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	305.83
Formulation	C15H16ClN3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
 SMILES: [Cl-].CC1=CC2=C(SC3=CC(C=CC3=N2)=[N+](C)C)C=C1N
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: Black Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Taha M, et al. Mater Sci Eng C Mater Biol Appl. 2013 Jul 1;33(5):2639-47.

Note: This product is for in vitro research use only