

Acedoben

Catalog No: #S0315


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Acedoben
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	179.17
Formulation	C ₉ H ₉ NO ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CC(=O)NC1=CC=C(C=C1)C(O)=O
Pathways: Others
Receptor: Others
Boiling pt: 439.6°C
Melting pt: 259-262°C
Solubility: Soluble in Water
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Franklin MR, et al. Xenobiotica. 1971 Mar;1(2):121-30.

Note: This product is for in vitro research use only