

## Citolone


Catalog No: #S0316

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Citolone
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	159.21
Formulation	C6H9NO2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
SMILES: CC(=O)NC1CCSC1=O  
Pathways: Others  
Receptor: Others  
Boiling pt: 427.3°C  
Melting pt: 109-111°C  
Solubility: DMSO: 10 mM  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. de Barrio M, et al. J Investig Allergol Clin Immunol. 1997 May-Jun;7(3):193-4.

Note: This product is for in vitro research use only