

Buflomedil hydrochloride

Catalog No: #S0317


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Buflomedil hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Adrenergic Receptor antagonist
Calculated MW	343.85
Formulation	C17H26ClNO4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System SMILES: Cl.COC1=CC(OC)=C(C(=O)CCCN2CCCC2)C(OC)=C1 Pathways: GPCR/G

Protein Receptor: α -adrenergic receptor Boiling pt: 454.5°C Melting pt: Solubility: Ethanol: 7 mg/mL (20.35 mM); Water: 69 mg/mL (200.66 mM); DMSO: 31 mg/mL (90.15 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Redondo P, et al. Exp Dermatol, 1997, 6(4), 186-194.

Note: This product is for in vitro research use only