

Ataractan

Catalog No: #S0318


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ataractan
Brief Description	Inhibitors
Purification	98.00%
Target Name	5-HT Receptor antagonist
Calculated MW	267.36
Formulation	C18H21NO
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:OC(C1CCNCC1)(C1=CC=CC=C1)C1=CC=CC=C1Pathways:GPCR/G Protein;

NeuroscienceReceptor:5-HTBoling pt:Melting pt:160-163CSolubility:Ethanol: 53 mg/mL (198.22 mM); Water: 2 mg/mL (7.48 mM); DMSO: 53 mg/mL (198.22 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Raeissi SD, et al. Pharm Res, 1997, 14(8), 1019-1025.

Note: This product is for in vitro research use only