

Ataractan

Catalog No: #S0318


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ataractan
Brief Description	Inhibitors
Purification	98.00%
Target Name	5-HT Receptor antagonist
Calculated MW	267.36
Formulation	C18H21NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system SMILES: OC(C1CCNCC1)(C1=CC=CC=C1)C1=CC=CC=C1 Pathways: GPCR/G Protein;

Neuroscience Receptor: 5-HT Boling pt: Melting pt: 160-163°C Solubility: Ethanol: 53 mg/mL (198.22 mM); Water: 2 mg/mL (7.48 mM); DMSO: 53 mg/mL (198.22 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Raeissi SD, et al. Pharm Res, 1997, 14(8), 1019-1025.

Note: This product is for in vitro research use only