

Indoprofen

Catalog No: #S0321


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Indoprofen
Brief Description	Inhibitors
Purification	98.00%
Target Name	COX inhibitor
Calculated MW	281.31
Formulation	C17H15NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CC(C(O)=O)C1=CC=C(C=C1)N1CC2=CC=CC=C2C1=OPathways:Immunology/Inflammation;
NeuroscienceReceptor:COXBoling pt:511.3C?Melting pt:Solubility:DMSO: 38mg/mLAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lunn MR, et al. Chem Biol. 2004 Nov;11(11):1489-93.

Note: This product is for in vitro research use only