

Fipexide

Catalog No: #S0322


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fipexide
Brief Description	Inhibitors
Purification	98.00%
Target Name	Dopamine Receptor agonist
Calculated MW	388.85
Formulation	C ₂₀ H ₂₁ CIN ₂ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: ClC1=CC=C(OCC(=O)N2CCN(CC3=CC4=C(OCO4)C=C3)CC2)C=C1
 Pathways: GPCR/G Protein;
 Neuroscience Receptor: Dopamine
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Marino A, et al. Pharmacol Res. 1990 Mar-Apr;22(2):179-87.

Note: This product is for in vitro research use only