

Isoxsuprine hydrochloride

Catalog No: #S0324


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Isoxsuprine hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Adrenergic Receptor agonist
Calculated MW	337.84
Formulation	C18H24ClNO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System SMILES: Cl.CC(COC1=CC=CC=C1)NC(C)C(O)C1=CC=C(O)C=C1 Pathways: GPCR/G

Protein Receptor: β -adrenergic receptor Boiling pt: 484.2°C Melting pt: Solubility: Soluble in Water Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Gozo EG Jr, Yebes RB. Chest. 1984 Nov;86(5):736-40.

Note: This product is for in vitro research use only