

Diperodon hydrochloride

Catalog No: #S0325


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Diperodon hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	433.93
Formulation	C ₂₂ H ₂₈ CIN ₃ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Other
 SMILES:O=C(NC1=CC=CC=C1)OCC(CN1CCCCC1)OC(=O)NC1=CC=CC=C1
 Pathways:Others
 Receptor:Others
 Boiling pt:
 Melting pt:502C
 Solubility:Soluble in DMSO
 Appearance:White Powder
 Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
 Stock solution can be stored below -20 °C for several month.

References

1.Calnan CD. Contact Dermatitis. 1980 Aug;6(5):367.

Note: This product is for in vitro research use only