

Diperodon hydrochloride

Catalog No: #S0325


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Diperodon hydrochloride |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | Others |
| Calculated MW | 433.93 |
| Formulation | C ₂₂ H ₂₈ CIN ₃ O ₄ |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Other
SMILES: O=C(NC1=CC=CC=C1)OCC(CN1CCCCC1)OC(=O)NC1=CC=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: 502°C
Solubility: Soluble in DMSO
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Calnan CD. Contact Dermatitis. 1980 Aug;6(5):367.

Note: This product is for in vitro research use only