

Fendiline hydrochloride

Catalog No: #S0327


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fendiline hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Calcium Channel inhibitor
Calculated MW	351.91
Formulation	C ₂₃ H ₂₆ ClN
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: Cl.CC(NCCC(C1=CC=CC=C1)C1=CC=CC=C1)C1=CC=CC=C1
 Pathways: Membrane transporter/Ion channel
 Receptor: Calcium Channel
 Boiling pt: 452.1°C
 Melting pt: Solubility: Water: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Sprgel W, et al. *Arzneimittelforschung*. 1977;27(3):571-4.

Note: This product is for in vitro research use only