

## Hexetidine

Catalog No: #S0328


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Hexetidine
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	339.61
Formulation	C21H45N3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:ImmunologySMILES:CCCCC(CC)CN1CN(CC(CC)CCCC)CC(C)(N)C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Water: 2.34 mg/LAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Afennich F, et al. Int J Dent Hyg. 2011 Aug;9(3):182-90.

Note: This product is for in vitro research use only