

Thiodiglycol


Catalog No: #S0329

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thiodiglycol
Brief Description	Inhibitors
Purification	97.00%
Target Name	AChR agonist
Calculated MW	122.19
Formulation	C4H10O2S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:CancerSMILES:OCCSCCOPathways:NeuroscienceReceptor:ADHBoiling pt:164-166CMelting pt:-16CSolubility:Water: 1 g/mL (at 20 C)Appearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Angerhofer RA, et al. Int J Toxicol. 2014 Sep-Oct;33(5):393-402.

Note: This product is for in vitro research use only