

Rosiglitazone

Catalog No: #S0334


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rosiglitazone
Brief Description	Inhibitors
Purification	98.00%
Target Name	PPAR agonist
Calculated MW	357.43
Formulation	C18H19N3O3S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism

systemSMILES:CN(CCOC1=CC=C(CC2SC(=O)NC2=O)C=C1)C1=CC=CC=N1Pathways:MetabolismReceptor:PPARBoiling pt:585CMelting

pt:153-155 CSolubility:DMSO: 71 mg/mL (198.64 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Pftzner A, et al. Vasc Health Risk Manag. 2007; 3(2):211-20.

Note: This product is for in vitro research use only