

Benzalkonium chloride

Catalog No: #S0336


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Benzalkonium chloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	N/A
Formulation	N/A
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: [Cl-].C[N+](C)(C)CC1=CC=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: >100°C
Melting pt: Solubility: Water: 10 mM
Appearance: Liquid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Baudouin C, et al. Prog Retin Eye Res. 2010 Jul;29(4):312-34.

Note: This product is for in vitro research use only