

## Bisotrizole

Catalog No: #S0340


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Bisotrizole
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	658.37
Formulation	C41H50N6O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: CC(C)(C)CC(C)(C)C1=CC(N2N=C3C=CC=C3=N2)=C(O)C(CC2=C(O)C(=CC(=C2)C(C)(C)CC(C)(C)C)N2N=C3C=CC=C3=N2)=C1 Pathways: Others Receptor: Others Boiling pt: 771.6°C Melting pt: 197-199°C Solubility: Chloroform: 35mg/mL Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Herzog B, et al. Int J Cosmet Sci. 2002 Jun;24(3):170-85.

Note: This product is for in vitro research use only