

## Miconazole nitrate

Catalog No: #S0344L


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Miconazole nitrate
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	479.15
Formulation	C <sub>18</sub> H <sub>15</sub> Cl <sub>4</sub> N <sub>3</sub> O <sub>4</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: O[N+](=O)([O-])=O.C1=CC(Cl)=C(COC(CN2C=CN=C2)C2=C(Cl)C=C(Cl)C=C2)C=C1 Pathways: Others Receptor: Others Boiling pt: 555.1 Melting pt: 170-185 Solubility: Ethanol: 1 mg/mL (2.08 mM); DMSO: 96 mg/mL (200.35 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Birsan M, et al. Rev Med Chir Soc Med Nat Iasi. 2014 Jan-Mar;118(1):225-31.

Note: This product is for in vitro research use only