

Estropipate


Catalog No: #S0345

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Estropipate
Brief Description	Inhibitors
Purification	98.86%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	436.56
Formulation	C22H32N2O5S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer SMILES: C[C@]12CC[C@H]3[C@H]([C@@H]1CCC2=O)CCC4=C3C=CC(=C4)OS(=O)(=O)O.C1CNCCN1 Pathways: Endocrinology/Hormones Receptor: ERBoling pt: Melting pt: Solubility: DMSO: 35mg/mL Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Brama M, et al. Mol Cell Endocrinol. 2007 Jan 29;264(1-2):102-8.

Note: This product is for in vitro research use only