

Phenol Red

Catalog No: #S0347


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Phenol Red
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	354.38
Formulation	C19H14O5S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism

systemSMILES:OC1=CC=C(C=C1)C1(OS(=O)(=O)C2=CC=CC=C2)C1=CC=C(O)C=C1Pathways:OthersReceptor:OthersBoiling pt:>300CMelting pt:300CSolubility:DMSO: 10 mMAppearance:Red PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Tamura, Z. Yakugaku Zasshi (in Japanese), 1997, 117 (10C11): 764C770.

Note: This product is for in vitro research use only