

Osalmid


Catalog No: #S0353

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Osalmid
Brief Description	Inhibitors
Purification	98.00%
Target Name	DNA/RNA Synthesis inhibitor
Calculated MW	229.23
Formulation	C13H11NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:OC1=CC=C(NC(=O)C2=CC=CC=C2O)C=C1Pathways:DNA Damage/DNA RepairReceptor:RRBoling
pt:350.8CMelting pt:179CSolubility:Soluble in Ethanol, DMSOAppearance:Black PowderRemark:For obtaining a higher solubility , please warm the
tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Sadanaga T, et al. J Chromatogr. 1981 Apr 10;223(1):243-6.

Note: This product is for in vitro research use only