

## Canrenone

Catalog No: #S0354


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Canrenone
Brief Description	Inhibitors
Purification	98.00%
Target Name	Glucocorticoid Receptor antagonist
Calculated MW	340.46
Formulation	C22H28O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Cardiovascular

SystemSMILES:C[C@]12CCC(=O)C=C1C=C[C@@H]3[C@@H]2CC[C@]4([C@H]3CC[C@@]45CCC(=O)O5)CPathways:Endocrinology/HormonesReceptor:Mineralocorticoid ReceptorBoiling pt:Melting pt:158-160CSolubility:DMSO: > 10 mMAppearance:Off-white Color PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Derosa G, et al. Ann Med. 2015 Feb;47(1):47-52.

Note: This product is for in vitro research use only