

Valpromide


Catalog No: #S0355

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Valpromide
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	143.23
Formulation	C ₈ H ₁₇ NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CCCC(CCC)C(N)=O
Pathways: Others
Receptor: 2-epoxide hydrolase; Limonene-1
Boiling pt: 274.6°C
Melting pt: 123-125°C
Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Lemoine P, et al. Eur Psychiatry. 2000 Nov;15(7):424-32.

Note: This product is for in vitro research use only