

Nitromide

Catalog No: #S0356


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nitromide
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	211.14
Formulation	C7H5N3O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: NC(=O)C1=CC(=CC(=C1)[N+](=O)[O-])[N+](=O)[O-]
 Pathways: Others
 Receptor: Others
 Boiling pt: 314.5°C
 Melting pt: 183-185°C
 Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Chubis AI, Terebinov GI. Veterinaria. 1972 Dec;12:73-4.

Note: This product is for in vitro research use only