

Potassium iodide

Catalog No: #S0357


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Potassium iodide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Tyrosinase
Calculated MW	166.01
Formulation	IK
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine systemSMILES:[K+].[I-]Pathways:Proteases/ProteasomeReceptor:TyrosinaseBoiling pt:184CMelting pt:113CSolubility:Water: 33 mg/mL (198.79 mM); DMSO: 33 mg/mL (198.79 mM)Appearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Pisarev MA, Aiello LO. Acta Endocrinol (Copenh). 1976 Jun;82(2):298-305.

Note: This product is for in vitro research use only