

Aminoguanidine hydrochloride


Catalog No: #S0358

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aminoguanidine hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	NOS inhibitor
Calculated MW	110.55
Formulation	CH6N4HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: Cl.NN=C(N)N
 Pathways: Immunology/Inflammation
 Receptor: iNOS
 Boiling pt: 261.4°C
 Melting pt: 162-166°C
 Solubility: DMSO: 10 mM
 Appearance: White Crystal
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Imming P, et al. Nat Rev Drug Discov. 2006 Oct; 5(10):821-34.

Note: This product is for in vitro research use only