

2-Aminoheptane


Catalog No: #S0359

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	2-Aminoheptane
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	115.22
Formulation	C7H17N
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CCCCC(C)NPathways:OthersReceptor:OthersBoiling pt:Melting pt:142-144CSolubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. COMER JP, KENNEDY EE. J Am Pharm Assoc Am Pharm Assoc (Baltim). 1956 Jul;45(7):454-5.

Note: This product is for in vitro research use only