

chlophenadione

Catalog No: #S0360


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	chlophenadione
Brief Description	Inhibitors
Purification	98.00%
Target Name	Vitamin inhibitor
Calculated MW	256.69
Formulation	C15H9ClO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
SMILES: C1C=CC=C(C=C1)C1C(=O)C2=CC=CC=C2C1=O
Pathways: Metabolism
Receptor: VKB
olig pt: 434.7
Melting pt: 142-144°C
Solubility: Water: 17.9 mg/L
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. King DA. Postgrad Med J. 1965 Dec;41(482):783-5.

Note: This product is for in vitro research use only