

## chlophenadione

Catalog No: #S0360


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	chlophenadione
Brief Description	Inhibitors
Purification	98.00%
Target Name	Vitamin inhibitor
Calculated MW	256.69
Formulation	C15H9ClO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:Metabolism systemSMILES:ClC1=CC=C(C=C1)C1C(=O)C2=CC=CC=C2C1=OPathways:MetabolismReceptor:VKBoiling pt:434.7°CMelting pt:142-144°CSolubility:Water: 17.9 mg/LAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. King DA. Postgrad Med J. 1965 Dec;41(482):783-5.

Note: This product is for in vitro research use only