

Morantel tartrate

Catalog No: #S0363


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Morantel tartrate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	370.42
Formulation	C18H24N2O7S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:CC1=C(SC=C1)/C=C/C2=NCCCN2C.C(C(C(=O)O)O)(C(=O)O)OPathways:OthersReceptor:OthersBoiling pt:334.6CMelting pt:Solubility:Soluble in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lanusse CE, et al. J Vet Pharmacol Ther. 1992 Jun;15(2):117-23.

Note: This product is for in vitro research use only