

Morantel tartrate

Catalog No: #S0363


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Morantel tartrate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	370.42
Formulation	C18H24N2O7S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: CC1=C(SC=C1)/C=C/C2=NCCCN2C.C(C(C(=O)O)O)(C(=O)O)O
 Pathways: Others
 Receptor: Others
 Boiling pt: 334.6°C
 Melting pt: Solubility: Soluble in DMSO
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Lanusse CE, et al. J Vet Pharmacol Ther. 1992 Jun;15(2):117-23.

Note: This product is for in vitro research use only