

Cefadroxil

Catalog No: #S0366


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cefadroxil
Brief Description	Inhibitors
Purification	98.00%
Target Name	Antibacterial
Calculated MW	363.39
Formulation	C16H17N3O5S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC1=C(N2[C@@H]([C@@H](C2=O)NC(=O)[C@@H](C3=CC=C(C=C3)O)N)SC1)C(=O)O Pathways: Microbiology & virology

Receptor: PBPs Boiling pt: Melting pt: Solubility: Soluble in Water Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Williamson R, et al. Antimicrob Agents Chem Others. 1980 Oct;18(4):629-37.

Note: This product is for in vitro research use only