

## Fluocinolone acetonide

Catalog No: #S0367


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Fluocinolone acetonide
Brief Description	Inhibitors
Purification	99.42%
Target Name	Others
Calculated MW	452.49
Formulation	C24H30F2O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine

systemSMILES:C[C@]12C[C@@H]([C@]3([C@H]([C@@H]1C[C@@H]4[C@]2(OC(O4)(C)C)C(=O)CO)C[C@@H](C5=CC(=O)C=C[C@@]53C)F)F)

OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Ethanol: 12 mg/mL (26.51 mM); DMSO: 90 mg/mL (198.89 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only