

Cinepazide maleate

Catalog No: #S0368


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cinepazide maleate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	533.57
Formulation	C22H31N3O5C4H4O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous

systemSMILES:COC1=CC(=CC(=C1OC)OC)/C=C/C(=O)N2CCN(CC2)CC(=O)N3CCCC3.C(=C/C(=O)O)\C(=O)OPathways:OthersReceptor:OthersBo
 ling pt:637.8C?Melting pt:Solubility:Water: 107 mg/mL (200.53 mM); DMSO: 107 mg/mL (200.53 mM)Appearance:White PowderRemark:For obtaining
 a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several
 month.

References

1. Pourrias B, et al. Therapie. 1974 Jan-Feb;29(1):29-41.

Note: This product is for in vitro research use only