

Ethaverine hydrochloride

Catalog No: #S0369


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ethaverine hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	MAO inhibitor
Calculated MW	431.96
Formulation	C ₂₄ H ₃₀ ClNO ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous

system SMILES: Cl.CCOC1=C(OCC)C=C(CC2=NC=CC3=CC(OCC)=C(OCC)C=C23)C=C1 Pathways: Metabolism Receptor: MAO Boling

pt: 527.1°C Melting pt: 187°C Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Lee SS, et al. Biol Pharm Bull. 2001 Jul;24(7):838-40.

Note: This product is for in vitro research use only