

## Gabexate mesylate

Catalog No: #S0372


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Gabexate mesylate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Serine Protease inhibitor
Calculated MW	417.48
Formulation	C17H27N3O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: CS(O)(=O)=O.CCOC(=O)C1=CC=C(OC(=O)CCCCN=C(N)N)C=C1 Pathways: Proteases/Proteasome Receptor: Serine Protease  
 Boiling pt: 508.6°C Melting pt: Solubility: Ethanol: 83 mg/mL (198.81 mM); Water: 12 mg/mL (28.74 mM); DMSO: 83 mg/mL (198.81 mM)  
 Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Erba F, et al. Biochem Pharmacol. 2001 Feb 1;61(3):271-6.

Note: This product is for in vitro research use only