

Atropine


Catalog No: #S0375

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Atropine
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR antagonist
Calculated MW	289.37
Formulation	C17H23NO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CN1C2CCC1CC(C2)OC(=O)C(CO)C1=CC=CC=C1
Pathways: Neuroscience
Receptor: mAChR
Boiling pt: 429.8°C
Melting pt: 115-118°C
Solubility: soluble in Water
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Huang XP, et al. J Pharmacol Exp Ther. 1998 Sep; 286(3):1129-39.

Note: This product is for in vitro research use only