

## Eprodinate disodium


Catalog No: #S0382

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Eprodinate disodium
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	248.19
Formulation	C <sub>3</sub> H <sub>6</sub> Na <sub>2</sub> O <sub>6</sub> S <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: [Na+].[Na+].[O-]S(=O)(=O)CCCS([O-])(=O)=O  
 Pathways: Others  
 Receptor: SAAB  
 Melting pt: >300°C  
 Solubility: Soluble in Water  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Bayes M, et al. Methods Find Exp Clin Pharmacol. 2007 Apr;29(3):231-45.

Note: This product is for in vitro research use only