

Prochlorperazine Maleate

Catalog No: #S0385


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Prochlorperazine Maleate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Dopamine Receptor agonist
Calculated MW	606.09
Formulation	C28H32ClN3O8S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous

system SMILES: CN1CCN(CC1)CCCN2C3=C(SC4=CC=CC=C4)C=CC(=C3)Cl.C(=C\C(=O)O)\C(=O)O.C(=C\C(=O)O)\C(=O)O
 Pathways: GPCR/G Protein; Neuroscience Receptor: D2
 Boiling pt: Melting pt: 175°C
 Solubility: Water: 15 mg/L (at 24°C)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37°C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

References

1. Hamik A, Peroutka SJ. Y Chem Others Pharmacol. 1989;24(5):307-10.

Note: This product is for in vitro research use only