

Prochlorperazine Maleate

Catalog No: #S0385


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Prochlorperazine Maleate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Dopamine Receptor agonist
Calculated MW	606.09
Formulation	C28H32ClN3O8S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous

systemSMILES:CN1CCN(CC1)CCCN2C3=C(SC4=CC=CC=C4)C=CC(=C3)Cl.C(=C\C(=O)O)\C(=O)O.C(=C\C(=O)O)\C(=O)OPathways:GPCR/G Protein; NeuroscienceReceptor:D2Boiling pt:Melting pt:175CSolubility:Water: 15 mg/L (at 24 C)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Hamik A, Peroutka SJ. Y ChemOthers Pharmacol. 1989;24(5):307-10.

Note: This product is for in vitro research use only