

## Dabigatran etexilate

Catalog No: #S0389


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

|                   |  |
|-------------------|--|
| Product Name      | Dabigatran etexilate                           |
| Brief Description | Inhibitors                                     |
| Purification      | 98.00%   |
| Target Name       | Thrombin inhibitor                             |
| Calculated MW     | 627.73   |
| Formulation       | C34H41N7O5                                     |
| Storage           | 3 years -20°C powder;2 years -80°C in solvent; |

## Images



## Product Description

## Research

Area: OtherSMILES: CCCCCOC(=O)/N=C(/C1=CC=C(C=C1)NCC2=NC3=C(N2C)C=CC(=C3)C(=O)N(CCC(=O)OCC)C4=CC=CC=N4)\NPathways: OthersReceptor: ThrombinBoiling pt: Melting pt: Solubility: Ethanol: 12 mg/mL (19.11 mM); DMSO: 126 mg/mL (200.72 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Squizzato A, et al. Intern Emerg Med. 2009 Dec;4(6):479-84.

Note: This product is for in vitro research use only