

Alkomide

Catalog No: #S0390

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alkomide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	200.58
Formulation	C ₇ H ₅ ClN ₂ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: NC(=O)C1=C(Cl)C=C(C=C1)[N+](=O)[O-]
 Pathways: Others
 Receptor: Others
 Boiling pt: 335.8°C
 Melting pt: 171-173°C
 Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Baron RR, et al. Poult Sci. 1966 Mar;45(2):411-8.

Note: This product is for in vitro research use only